## **Torque Settings**




3.5 mm Torque to: 35 Ncm Replace Select 4.3 mm Torque to: 35 Ncm

5.0 mm Torque to: 35 Ncm 6.0 mm Torque to: 35 Ncm

Nobel Active 3.5 mm Torque to: 35 Ncm

4.3 mm Torque to: 35 Ncm 5.0 mm Torque to: 35 Ncm











Certain 3.4 mm Torque to: 20 Ncm 4.1 mm Torque to: 20 Ncm

5.0 mm Torque to: 20 Ncm

\*20 Ncm applies to gold AND titanium screws



External

3.4 mm Torque to: square screw 35 Ncm - Hexed screw 20 Ncm

4.1 mm Torque to: square screw 35 Ncm - Hexed screw 20 Ncm

5.0 mm Torque to: square screw 35 Ncm - Hexed screw 20 Ncm







3.5 mm Torque to: 35 Ncm Tissue Level 4.8 mm Torque to: 35 Ncm

4.8 mm Torque to: 35 Ncm 6.5 mm Torque to: 35 Ncm

3.3 mm Torque to: 35 Ncm


Bone Level 4.1 mm Torque to: 35 Ncm

4.6 mm Torque to: 35 Ncm



Tapered Screw Vent

3.5 mm Torque to: 30 Ncm 4.5 mm Torque to: 30 Ncm 5.7 mm Torque to: 30 Ncm





Camlog

3.3 mm Torque to: 20 Ncm 3.8 mm Torque to: 20 Ncm 4.3 mm Torque to: 20 Ncm 5.0 mm Torque to: 20 Ncm 6.0 mm Torque to: 20 Ncm





